

Data Sheet

IBA Headquarters IBA GmbH

Rudolf-Wissell-Str. 28 37079 Goettingen Germany Tel. +49 (0) 551-5 06 72-0 Fax +49 (0) 551-5 06 72-181

IBA US contact information Fax 1-888-531-6813

Version:

E-mail: info@iba-lifesciences.com http://www.iba-lifesciences.com

2.2 Revision Date: 04.03.2020

pESG-IBA23

Cat. No.: 5-4423-001

Lot No.: 4423-

Description Affinity tag	StarGate® Acceptor Vector designed for high-level stable and non-replicative transient expression in most mammalian hosts containing the following elements: • Human cytomegalovirus (CMV) immediate-early promoter for high-level expression in a wide range of mammalian cells • Neomycin resistance gene for selection of stable cell lines • Episomal replication in cell lines that are latently infected with SV40 or that express the SV40 large T antigen (e.g. COS-1, COS-7) • The expressed recombinant protein will contain two affinity tags:		
Affinity tag	 The recombinant protein will contain two affinity tags: Strep-Tactin affinity tag (Strep-tag® II) for the purification of recombinant protein via Strep-Tactin resins. The Strep-tagII is fused to the C-terminus of the recombinant protein. GST-tag (Glutathione-S-Transferase) for the purification of recombinant protein. The affinity tag is fused to the N-terminus of the recombinant protein. After purification the GST may be removed by digesting with PreScission™ Protease. 		
Resistance	Ampicillin		
Form	5 μg, dissolved in 20 μl TE buffer, pH 8,0: 10 mM Tris-HCl, 1 mM EDTA		
Concentration	250 ng/μl		
Stability	12 months after shipping		
Storage	recommended: 2-8 °C for frequent usage, -20 °C for long-term storage		
Shipping	room temperature		
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.		

Note: The sequences have been compiled from information in the sequence database, published literature, and other sources, together with partial sequences obtained by IBA, however, the vectors have not been completely sequenced.

PreScission™ Protease is a trademark of GE HEALTHCARE

For research use only


Important licensing information

This product is based on StarGate, Strep-tag, GST-tag and CMV promoter technologies covered by intellectual property (IP) rights and on completion of the sale IBA grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at http://www.iba-lifesciences.com/patents.html or upon inquiry at info@iba-lifesciences.com or at IBA GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.

Trademark information

The owners of trademarks marked by "®" or "TM" are identified at http://www.iba-lifesciences.com/patents.html. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Expression cassette of pESG-IBA23




Expression cassette of pESG-IBA23, continued

LacP/Z cassette = contains LacZ alpha fragment under control of a separate promoter, which allows alpha

complementation of LacZ mutations such as $LacZ\Delta M15$ as in *E. coli* DH5 α or TOP10.

your protein =

after StarGate cloning using Esp3l your gene of interest will be located here



Features	from bp	to bp	Sequencing primer
f1 origin	259	687	ESG-Primer-for (Cat. No. 5-0000-121)
SV40 ori	692	1035	
Neomycin resistance gene	1097	1891	5'- GAGAACCCACTGCTTACTGGC -3'
ColElori	2637	3222	
Ampicillin resistance gene	3393	4253	ESG-Primer-rev (Cat. No. 5-0000-122)
CMV promoter	4621	5208	E3G-Fillier-lev (Cat. No. 3-0000-122)
forward primer binding site	5221	5241	5'- TAGAAGGCACAGTCGAGG -3'
GST-tag	5284	5937	
PreScission™ protease site (PCS)	5938	5985	
LacZ alpha fragment	6214	6615	
Strep-tagll	6679	6711	
reverse primer binding site	6773	6790	
polyA signal sequence	1	213	
total vector length		6790	



Go digital and help the environment. Please download all up-to-date manuals, protocols and other material from http://www.iba-lifesciences.com.