



5-HTP (5-Hydroxytryptophan)

Data Sheet

| | | | |
|----------------------------|--|--------------------|---|
| Catalog Number: | RA24446 | Host: | Rabbit |
| Product Type: | Whole Serum | Species | Rat |
| Immunogen Sequence: | 5-HTP coupled to bovine serum albumin (BSA) with paraformaldehyde. | Reactivity: | |
| | | Format: | Lyophilized. 100 ul with <0.09% sodium azide as a preservative. |

Applications: Immunohistochemistry: 1:1,000-1:2,000 in PBS/0.3% Triton X-100 - Bn-AV/HR.

Optimal dilution will vary depending upon fixation, labeling technique and/or detection system; therefore, a dilution series is recommended. Staining is completely eliminated by pretreatment of the diluted antibody with 5 µg of LHRH per mL of diluted antiserum.

Reconstitution: Do not reconstitute until ready to use since the product is most stable when lyophilized. The product does not need to be kept cooled during shipping. For long-term storage, store lyophilized antibody until ready to use at -15° C or lower. Reconstitute with 100 µL of distilled or deionized water. If desired, dilute with PBS or Tris buffer at a dilution no higher than 1/10

Storage: After reconstitution, use immediately or refrigerate at 2°-8° C up to 2 days. For long-term storage, appropriately aliquot antibody to avoid repeated freeze/thaw cycles and freeze at -15° C or lower

Application Notes

Description/Data:

5-Hydroxytryptophan (5-HTP), also known as oxitriptan (INN), is a naturally occurring amino acid and chemical precursor as well as metabolic intermediate in the biosynthesis of the neurotransmitters serotonin and melatonin from tryptophan. 5-Hydroxytryptophan is decarboxylated to serotonin (5-hydroxytryptamine or 5-HT) by the enzyme aromatic-L-amino-acid decarboxylase with the help of vitamin B6.[4] This reaction occurs both in nervous tissue and in the liver. 5-HTP crosses the blood-brain barrier,[6] while 5-HT does not. Since 5-HTP increase levels of serotonin, it is used to treat depression, sleep disorders and pain.

Tissue Preparation:

10 µm cryostat or 50 µm vibratome.

- Fixation: 4% paraformaldehyde in 0.1M Phosphate buffer, pH 7.4; 500 mL over ~ 20 min.
- Post Fixation: 1.5 hour at 4° C in 4% paraformaldehyde in 0.1M phosphate buffer, pH 7.4.
- Note: If needed, low levels of glutaraldehyde (0.1-0.3%) may be used in conjugation with paraformaldehyde

Immunohistochemistry



The antibody has significant biotin-avidin/HRP staining at a 1/1,000 - 1/2,000 dilution in rat raphe nuclei. Optimal dilution will vary depending upon fixation, labeling technique and/or detection system; therefore, a dilution series is recommended. The specificity of the antiserum was evaluated using a model system of gelatin-indole plugs by a method similar to published procedures (Shipper and Tilders, 1983). Results showed that the 5-HTP antibody dose dependently stained 5-HTP but did not stain any concentration of 5-HT or 5-HIAA. The antiserum was also tested by pre-adsorption with indole/paraformaldehyde/BSA conjugates. Staining was completely blocked by pre-adsorption with 5-HTP conjugate and unaffected by 5-HIAA or 5-HT conjugate.

FOR RESEARCH USE ONLY

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.-V2/08/2012

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail: pshuster@neuromics.com



FOR RESEARCH USE ONLY

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.-V2/08/2012

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail: pshuster@neuromics.com