NEUROMICS

Calcineurin A

Data Sheet

Catalog Number:	MO15055	Host:	Mouse
Product Type:	Protein G purified rat IgG ₁ (302202).	Species Reactivity:	Human, Rat, Mouse
lmmunogen Sequence:	Hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a rat immunized with purified bovine brain Calcineurin holoenzyme containing both the A and B subunits.	Format:	Liquid 1mg/ml Solution in phosphate-buffered saline (PBS) with 5% Trehlose
Applications:	Western Blot: 1 µg/mL		
	Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.		
Storage:	Antibody can be aliquotted and stored frozen at -20° C to -70° C in a manual defrost freezer for size months without detectable loss of activity. The antibody can be stored at 2° - 8° C for 1 month without detectable loss of activity. <i>Avoid repeated freeze-thaw cycles.</i>		

Application Notes

Specificity:

The antibody detects endogenous human, mouse, and rat Calcineurin A at 59 kDa in Western blots.

Western Blot

Blotting BufferBlocking SolutionAntibody Solution25 mM Tris, pH 7.45% nonfat dry milk in Blotting Buffer5% nonfat dry milk in Blotting Buffer0.15 M NaClAdjust pH to 7.45% nonfat dry milk in Blotting Buffer0.1% Tween 20Adjust pH to 7.45% nonfat dry milk in Blotting Buffer

1. Transfer the electrophoresed proteins to Immobilon-P membrane (Millipore) and incubate the membrane for 1 hour at room temperature in Blocking Solution.

Incubate the membrane overnight at 4° C in Antibody Solution containing 1 µg/mL anti-human/mouse/rat Calcineurin A.
Wash the membrane at room temperature for 1 hour with 5 or more changes of Blotting Buffer. Changing the membrane containers often reduces background.

4. Incubate the membrane at room temperature for 1 hour in Antibody Solution containing a 1:1,000 dilution of HRP-

conjugated goat anti-rat IgG (Zymed).

5. Wash the membrane for 1 hour with 5 or more changes of Blotting Buffer.

6. Detect with WesternGlo Chemiluminescent Detection Substrate or equivalent.

Cell Lysates for Western Blottings: To prepare total cell lysates, cells are solubilized in hot 2x SDS gel sample buffer (20 mM dithiothreitol, 6% SDS, 0.25 M Tris, pH 6.8, 10% glycerol, 10 mM NaF, and bromophenyl blue) at 2×10^{7} - 1 x 10 cells per mL. The extracts are heated in a boiling water bath for 5 minutes and then sonicated with a probe sonicator with 3 - 4 bursts of 5 - 10 seconds each. Samples are diluted with 1x SDS sample buffer to the desired concentration.

Description/Data:

Calcineurin A, also known as PP2B and PPP3CA, is the 59 kDa catalytic subunit of the calcium/calmodulin-dependent protein phosphatase. When activated by calcium in the presence of the regulatory B subunit and calmodulin, Calcineurin A selectively removes phosphates from serine and threonine residues on target proteins. Although ubiquitously expressed, Calcineurin levels are highest in brain, where the phosphatase plays a role in the formation and retention of memories.

FOR RESEARCH USE ONLY

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RSKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.-V2/08/2012

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439 phone 866-350-1500 • fax 612-677-3976 • e-mail: <u>pshuster@neuromics.com</u>

Image: Recombinant human Calcineurin A (1 ng) and lysates prepared from human HT29, mouse TS1, and rat Nb2-11 cells were resolved by SDS-PAGE, transferred to an Immobilon-P membrane, and immunoblotted with 1 μ g/mL anti-Calcineurin A, as described on datasheet. A ten second exposure to film is shown.



FOR RESEARCH USE ONLY

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RSKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.-V2/08/2012

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439 phone 866-350-1500 • fax 612-677-3976 • e-mail: <u>pshuster@neuromics.com</u>