

## Data Sheet

# Strep-Tactin XT<sup>®</sup> 4Flow<sup>®</sup>

## High Capacity

### 50% suspension

Cat. No.: 2-5030-002, 2-5030-010,  
2-5030-025, 2-5030-500

Lot No.: 5030-

Version: 1.0  
 Revision Date: 02.03.2019

<b>Description</b>	Strep-Tactin <sup>®</sup> XT high capacity resin for the purification of Strep-tag <sup>®</sup> II and Twin-Strep-tag <sup>®</sup> fusion proteins. Strep-Tactin <sup>®</sup> XT is a streptavidin variant with optimized binding properties for Twin-Strep-tag <sup>®</sup> and Strep-tag <sup>®</sup> II fusion proteins.
<b>Support</b>	4% agarose, crosslinked
<b>Form</b>	50% suspension in buffer, pH 8.0: 100 mM Tris/HCl pH 8.0, 1 mM EDTA, 150 mM NaCl, 0.02% sodium azide
<b>Dynamic Binding Capacity</b>	16 mg protein/ml resin. Dynamic binding capacity was determined with 1 mg/ml mCherry-Twin-Strep-tag <sup>®</sup> (30 kDa) at a flow rate of 0.5 ml/min. Please note: Binding capacity is protein dependent.
<b>Stability</b>	6 months after shipping
<b>Storage</b>	recommended: 2–8 °C
<b>Shipping</b>	room temperature
<b>Hazards</b>	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.

<b>Application</b>	<b>Important note:</b> For efficient purification results we strongly recommend using column purification instead of batch applications. Even a pre-incubation of resin and lysate prior to filling the resin into a column will lead to decreased protein yields. Further, prolonged batch incubations increase the risk of proteolytic degradation of the target protein including cleavage of the tag. Batch purification should be performed using Twin-Strep-tag <sup>®</sup> in combination with MagStrep“type3”XT beads only.
<b>Elution</b>	Strep-Tactin <sup>®</sup> XT Elution Buffer with biotin (10 x Buffer BXT, Cat. No. 2-1042-025); 1 x Buffer BXT: 100 mM Tris-HCl pH 8.0, 150 mM NaCl, 1 mM EDTA, 50 mM biotin
<b>Regeneration</b>	Regenerate the column by using Strep-Tactin <sup>®</sup> XT Regeneration Buffer (3 M MgCl <sub>2</sub> , Cat. No. 2-1045-250)

\* Voss, S. & Skerra, A. (1997) Mutagenesis of a flexible loop in streptavidin leads to higher affinity for the Strep-tag II peptide and improved performance in recombinant protein purification. *Protein Eng.* 10, 975-982.

### For research use only

#### Important licensing information

This product is based on Strep-Tactin XT technology covered by intellectual property (IP) rights and on completion of the sale IBA grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at [info@iba-lifesciences.com](mailto:info@iba-lifesciences.com) or at IBA GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.

#### Trademark information

The owners of trademarks marked by “®” or “TM” are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.



Please download all up-to-date manuals, protocols and other material from <https://www.iba-lifesciences.com>.