



**IBA Lifesciences GmbH**  
Rudolf-Wissell-Str. 28  
37079 Goettingen  
Germany  
Tel.: +49 (0) 551-5 06 72-0  
E-mail: [info@iba-lifesciences.com](mailto:info@iba-lifesciences.com)  
[www.iba-lifesciences.com](http://www.iba-lifesciences.com)

## Data Sheet

# GFP-Twin-Strep-tag control protein

Cat. No.: 2-1007-105

Version: 4.1  
Revision Date: 04.05.2023

<b>Description</b>	<b>GFP-Twin-Strep-tag is a control protein consisting of the green fluorescent protein (GFP) C-terminally tagged with the Twin-Strep-tag.</b>
<b>MW</b>	30 kDa
<b>Form</b>	Dissolved in 1x Buffer W
<b>Concentration</b>	1 mg/ml
<b>Amount</b>	500 µl
<b>Stability</b>	12 months after shipping
<b>Storage</b>	-25 °C to -15 °C
<b>Shipping</b>	Room temperature
<b>Hazards</b>	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.
<b>Application</b>	It is suitable for test purifications, detections, and immobilizations with Strep-Tactin® or Strep-Tactin®XT.

### For research use only

#### Trademark information

The owners of trademarks marked by “®” or “TM” are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

#### Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale IBA Lifesciences grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at [info@iba-lifesciences.com](mailto:info@iba-lifesciences.com) or at IBA Lifesciences GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.