

NEUROMICS



Human IPSC Astrocyte Growth Medium

Catalog #: PGB005

Size: 500 ml

Storage: Store the medium at 2-8°C or -20°C for long term

Shelf Life: 6 Months

Sterilization: 0.2 µm sterile filtered

Format: Liquid

GENERAL INFORMATION

Human IPSC Astrocyte Growth Medium is designed for in vitro culture of iPSC-derived astrocytes. This medium is fully supplemented. This medium is sterile-filtered and is ready for use.

The medium has been tested and is negative for bacteria, yeast, fungi, and mycoplasma.

Product is for Research use only. Our products are not authorized for human use, in vitro diagnostic procedures, or for therapeutic procedures.

Storage Condition: Media should be stored at 2 to -20°C. If taking aliquots from the original bottle is necessary, dispense the contents and re-freeze in small aliquots to avoid repeated freezing-thawing cycles.

FOR RESEARCH USE ONLY

NEUROMICS REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. v1-09809

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail pshuster@neuromics.com