



Catalog Number:	CH23003	Host:	Chicken
Product Type:	Affinity Purified	Species	Human, Mouse
Immunogen Sequence:	Peptide corresponded to a region of the NSE gene product shared between the rat (AAA41119, NCBI) and human (NP_001966, NCBI) gene products.	Reactivity:	
		Format:	Liquid PBS (pH 7.2; 10 mM; isotonic 0.9%, w/v) with sodium azide (0.02%, w/v).
Applications:	Immunocytochemistry: 1:100 -1:200 Immunohistochemistry: 1:100 -1:200 Western blot: 1:250 -1:500		
	Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.		
Storage:	Store at 4°C in the dark. Under these conditions, the antibody should have a shelf life of at least 12 months (provided they remain sterile). Do not freeze the antibody unless you want to store it for longer periods of time. Note, however, that each time an antibody preparation is frozen, about half its binding activity is lost.		

Application Notes