

NEUROMICS



TRPM8

Data Sheet

Catalog Number:	RA19043	Host:	Rabbit
Product Type:	Affinity Purified	Species	Human
Immunogen Sequence:	Reacts with residues 278-292 RNQLEKYISERTIQD and the C terminus sequence NDLKGLLKEIANKIK of the human TRPM8 Protein.	Reactivity:	
Applications:	Western Blot: 1:500-1000	Format:	100 µg in 100 µl (1 mg/ml) in PBS containing 0.02% sodium azide.
Storage:	Dilutions listed as a recommendation. Optimal dilution should be determined by investigator. Maintain at +2-8°C for 3 months or at -20°C for longer periods. Stable for 1 year. <i>Avoid repeated freeze-thaw cycles.</i>		

Application Notes

FOR RESEARCH USE ONLY

NEUROMICS REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. v3-01/13

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail pshuster@neuromics.com