



---

<b>Catalog Number:</b>	MO16101	<b>Host:</b>	Mouse
<b>Product Type:</b>	Protein G purified IgG <sub>1</sub> <b>Clone:</b> 4F7.2	<b>Species Reactivity:</b>	Human
<b>Immunogen Sequence:</b>	Recombinant human LIF	<b>Format:</b>	Liquid 0.02M PBS, 0.25M NaCl pH 7.6 with 0.1% Sodium Azide
<b>Applications:</b>	<b>Western blot:</b> 1 µg/mL Dilutions listed as recommendation. Optimal dilution should be determined by the end user.		
<b>Storage:</b>	Maintain refrigerated at 2-8°C in undiluted aliquots for up to 12 months.		

---

### Application Notes

This antibody is reactive to recombinant human Leukemia Inhibitory Factor (LIF) and to endogenous protein. Cross reactivity was observed against recombinant murine LIF. No reactivity was observed with recombinant rat LIF