



Catalog Number:	RA24506	Host:	Rabbit
Product Type:	Affinity purified	Species Reactivity:	Rat, Human, Mouse
Immunogen Sequence:	Synthetic peptide sequence corresponding to amino acids 356-382 of the rat NPY Y1 receptor.	Format:	Liquid 100 ul ≤ 0.09% sodium azide
Applications:	Immunohistochemistry: 1:500-1:1000 Immunocytochemistry: 1:500 Western Blot: 1:500-1:1000		The antiserum is provided as 100 µL of affinity purified liquid containing 1% BSA. Reconstitution is not required. Concentration: 300 ug/ml
Storage:	Dilutions listed as a recommendation. Optimal dilution should be determined by investigator. Maintain at +2-8°C for 3 months or at -20°C for longer periods. Stable for 1 year. <i>Avoid repeated freeze-thaw cycles.</i>		

Application Notes

To improve staining intensity it may be necessary to employ HRP-DAB (Vector Elite ABC kit) detection or enhance fluorescence signal using Tyramide Signal Amplification procedure.