

SHANK1a

Data Sheet

Catalog Number: RA19015

Product Type: Rabbit IgG

Immunogen Sequence: ERLREY PQSFEKGV of the N terminus of SHANK1a

Format: Affinity Purified lyophilized

Species Reactivity: Human, Mouse, Rat

Applications: Western blotting

Dilutions: WB 1:500 -1:1000

Storage: Maintain at +2-8 degrees C for 3 months or at -20 degrees C for longer periods. Stable for 1 year. Avoid repeated freeze-thaw cycles.

Application Note

Anti-SHANK1a antibody is at a concentration of 1 mg/ml in PBS containing 0.02% sodium azide.