



Catalog Number: RA19009

Host: Rabbit

Product Type: Affinity Purified Antibody

Species Reactivity: Human

Immunogen Sequence: GPGSRARAAGARGC

Format: Liquid. 100 ug in 100 ul (1 mg/ml) in PBS containing 0.02% sodium azide


Applications: Western blot: 1:500

Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.

Storage: Maintain at +2-8°C for 3 months or at -20°C for longer periods. Stable for 1 year. *Avoid repeated freeze-thaw cycles.*

Publications: [Deborah H. Damon, Jaclyn A. teRiele, and Stephen B. Marko. Vascular-derived artemin: a determinant of vascular sympathetic innervation? Am J Physiol Heart Circ Physiol 293: H266-H273, 2007. First published March 2, 2007; doi:10.1152/ajpheart.00859.2006 0363-6135/07.](#)

Application Notes



Western blot analysis of human Artemin in dorsal root ganglia / CNS preparation with anti-Artemin at a dilution of 1:500.

FOR RESEARCH USE ONLY

NEUROMICS REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. v3-01/13

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail pshuster@neuromics.com