



Catalog Number:	GT15105	Host:	Goat
Product Type:	Affinity purified	Species Reactivity:	Rat
Immunogen Sequence:	Purified, <i>E. coli</i> -derived, recombinant rat interleukin 1 beta (rrIL-1 β)	Format:	Liquid 1mg/ml Solution in phosphate-buffered saline (PBS) with 5% Trehlose
Applications:	Immunohistochemistry : 0.5-5 μ g/mL Western Blot: 0.1 - 0.2 μ g/mL ELISA Capture: 0.2 – 0.8 μ g/mL		

Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.

Storage: Antibody can be aliquotted and stored frozen at -20° C to -70° C in a manual defrost freezer for six months without detectable loss of activity. The antibody can be stored at 2° - 8° C for 1 month without detectable loss of activity. *Avoid repeated freeze-thaw cycles.*

Application Notes

Specificity

This antibody has been selected for its ability to neutralize the biological activity of rrIL-1 β and for use as a capture antibody in rat IL-1 β sandwich ELISAs. It will also neutralize the biological activity of rmlL-1 β , but will not neutralize the biological activity of rhIL-1 β , rhIL-1 α or rmlL-1 α . It can also be utilized for immunohistochemistry and western blotting.

Western blot

This antibody can be used at 0.1 – 0.2 μ g/mL with the appropriate secondary reagents to detect rat IL-1 β . The detection limit for rrIL-1 β is approximately 5 ng/lane and 2 ng/lane under non-reducing and reducing conditions, respectively. In Western blots, this antibody shows approximately 50% cross-reactivity with rmlL-1 β , rhIL-1 β and rpIL-1 β .

Immunohistochemistry

This antibody can be used at 0.5 - 5 μ g/mL with the appropriate secondary reagents to detect rat IL-1 β in cultured cells or tissue sections.

ELISA capture

This product can be used as a capture reagent in a rat IL-1 β sandwich immunoassay in combination with biotinylated rat IL-1 β detection antibody and recombinant rat IL-1 β as the standard. The suggested coating concentration range is 0.2 - 0.8 μ g/mL and should be titrated to determine the optimal concentration. In this format, less than 1% cross-reactivity with rmlL-1 β and less than 0.5% cross-reactivity with rhIL-1 β and rpIL-1 β is observed.

FOR RESEARCH USE ONLY

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RSKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

03/09v1

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail pshuster@neuromics1.com

Neutralization of Rat IL-1 β bioactivity - The exact concentration of antibody required to neutralize rrIL-1 β activity is dependent on the cytokine concentration, cell type, growth conditions and the type of activity studied. To provide a guideline, R&D Systems has determined the neutralization dose for this antibody under a specific set of conditions. The **Neutralization Dose₅₀ (ND₅₀)** for this antibody is defined as that concentration of antibody required to yield one-half maximal inhibition of the cytokine activity on a responsive cell line, when that cytokine is present at a concentration just high enough to elicit a maximum response.




Figure 1: Rat IL-1 β stimulates ³H-thymidine incorporation by murine T-helper D10.G4.1 cells in a dose-dependent manner (Symons, J.A. *et al.*, 1987, in *Lymphokines and Interferons, A Practical Approach*, IRL Press, M.J. Clemens, A.G. Morris and A.J.H. Gearing, eds. p. 272). The ED₅₀ for this effect is typically 1 - 3 ng/mL.

Figure 1




Figure 2: Approximately 0.5 - 2.0 μ g/mL of the antibody will neutralize 50% of the bioactivity due to 10 ng/mL of rat IL-1 β . The ND₅₀ for this lot of anti-rat IL-1 β antibody was determined to be approximately 0.5 - 2 μ g/mL in the presence of 10 ng/mL of rrIL-1 β , using the murine T-helper cell line, D10.G4.1. The specific conditions are described in the figure legends.

Figure 2

FOR RESEARCH USE ONLY

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

03/09v1

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail pshuster@neuromics1.com