

# NEUROMICS



## Neuropeptide Q/NPQ

## Data Sheet

---

|                            |   |                    |   |
|----------------------------|---|--------------------|---|
| <b>Catalog Number:</b>     | RA14131   | <b>Host:</b>       | Rabbit  |
| <b>Product Type:</b>       | Whole Serum   | <b>Species</b>     | Rat, Human                                    |
| <b>Immunogen Sequence:</b> | NWTPQAMLYLKGAQ-CONH2 Synthetic Peptide from the C-Terminal of NPQ.  | <b>Reactivity:</b> |   |
|                            |   | <b>Format:</b>     | 100 ul. in PBS containing 0.05% sodium azide. |
| <b>Applications:</b>       | Immunohistochemistry: 1:100-600<br>Dilutions listed as a recommendation. Optimal dilution should be determined by investigator. |                    |   |
| <b>Storage:</b>            | Maintain at +2-8°C for 3 months or at -20°C for longer periods. Stable for 1 year. <i>Avoid repeated freeze-thaw cycles.</i>    |                    |   |

---

### Application Notes

To improve staining intensity it may be necessary to employ HRP-DAB (Vector Elite ABC kit) detection or enhance fluorescence signal using Tyramide Signal Amplification procedure.

---

FOR RESEARCH USE ONLY

Neuromics' reagents are for in vitro and certain non-human in vivo experimental use only and not intended for use in any human clinical investigation, diagnosis, prognosis, or treatment. We disclaim all liability in connection with the use of the information contained herein or otherwise, and all such risks are assumed by the user.

v1