

Neuropeptide FF2 Receptor

Data Sheet

Catalog Number: RA19018

Product Type: Rabbit IgG

Immunogen Sequence: N terminus (MNEKWDTNSSSENWHPI) and C terminus (ELVMEELKETTNSSEI) of human NPFF2 protein.

Format: Affinity Purified liquid

Species Reactivity: Human

Applications: Western blotting

Dilutions: WB 1:500-1:1000

Storage: Maintain at +2-8 C for 3 months or at -20 C for longer periods. Stable for 1 year. Avoid repeated freeze-thaw cycles