



---

<b>Catalog Number:</b>	CH23005	<b>Host:</b>	Chicken
<b>Product Type:</b>	Affinity Purified	<b>Species Reactivity:</b>	Human, Mouse
<b>Immunogen Sequence:</b>	Peptide corresponded to a region of the Tuj 1 gene product shared between the rat (AAM28438, NCBI) and human (AAL28094, NCBI) gene products.	<b>Format:</b>	Liquid PBS (pH 7.2; 10 mM; isotonic 0.9%, w/v) with sodium azide (0.02%, w/v).
<b>Applications:</b>	<b>Immunocytochemistry:</b> 1:100 -1:200 <b>Immunohistochemistry:</b> 1:100 -1:200 <b>Western blot:</b> 1:250 -1:500 Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.		
<b>Storage:</b>	<b>Store at 4°C in the dark.</b> Under these conditions, the antibody should have a shelf life of at least 12 months (provided they remain sterile). Do not freeze the antibody unless you want to store it for longer periods of time. Note, however, that each time an antibody preparation is frozen, about half its binding activity is lost.		

---

### Application Notes