



**Catalog Number:** RA19023

**Host:** Rabbit

**Product Type:** Affinity Purified Antibody

**Species Reactivity:** Human; Mouse; Rat and Primate

**Immunogen Sequence:** Reacts with the residues 780-796 [YDDYEDDPKSDYSSENE] of human and mouse SOX6. Sequence 100% conserved between human, mouse, rat and chimpanzee.

**Format:** Liquid. 100 ug in 100 ul (1 mg/ml) in PBS containing in phosphate-buffered saline (PBS) with 5% Trehlose.

**Applications:** Immunohistochemistry: 1:100-1:300  
Western Blot: 1:1,000-1:2,000  
Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.

**Storage:** Antibody can be aliquotted and stored frozen at -20° C to -70° C in a manual defrost freezer for sixmonths without detectable loss of activity. The antibody can be stored at 2° - 8° C for 1 month without detectable loss of activity. *Avoid repeated freeze-thaw cycles*

### Application Notes