



Catalog Number: RA19000

Host: Rabbit

Product Type: Affinity purified antibody

Species Reactivity: Human

Immunogen Sequence: LLLLTSEIDLPKRRA
a.a (9-24 of the Humanin protein)

Format: **Liquid**
1mg/ml in PBS containing
0.02% sodium azide.

Applications: **Western blotting** 1:500
Immunoprecipitation Not Tested
Immunocytochemistry Not Tested

Dilutions listed only as a recommendation. Optimal dilution should be determined by investigator.

Storage: Maintain at +2-8°C for 1 month or at -20°C for longer periods. Liquid may have become trapped in top of vial during shipping. Centrifugation of vial is recommended before opening. Stable for at least 6 months at -20°C. Repeated freeze/thaw cycles compromise the integrity of the antibody.

Application Notes

FOR RESEARCH USE ONLY

NEUROMICS REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

www.neuromics.com

Neuromics Antibodies • 11200 Hampshire Ave South • Bloomington, MN 55438
phone 507-645-8020 • fax 612-677-3976 • e-mail pshuster@neuromics1.com