

SHANK1a C-terminus

Data Sheet

Catalog Number: RA19016

Product Type: Rabbit IgG

Immunogen Sequence: SGPIYPGLFDIRSS of the C-terminus of SHANK1a

Format: Affinity Purified liquid

Species Reactivity: Rat

Applications: Western blotting

Dilutions: WB 1:500 - 1:1000

Storage: Maintain at +2-8 degrees C for 3 months or at -20 degrees C for longer periods. Stable for 1 year. Avoid repeated freeze-thaw cycles.

Application Note

Anti- SHANK1a C-terminus antibody is at a concentration of 1 mg/ml in PBS containing 0.02% sodium azide.