



Catalog Number: RA25029

Host: Rabbit

Product Type: Affinity Purified Antibody

Species Reactivity: Rat, Human, Mouse

Immunogen Sequence: Recombinant rat protein expressed in bacteria.

Format: 100 ul (1 mg/ml) in PBS containing 0.02% sodium azide and 15% glycerol.

Applications: Western Blot : 1:1000

Dilutions listed only as a recommendation. Optimal dilution should be determined by investigator.

Storage: Store frozen. Aliquot as undiluted antisera and immediately place at -20°C. Antisera may have become trapped in top of vial during shipping. Centrifugation of vial is recommended before opening. Stable for at least 6 months at -20°C. Repeated freeze/thaw cycles compromise the integrity of the antiserum.


Application Notes

Recommended for Western Blot Applications. A band is seen at ~100 kDa. The observed molecular weight corresponds to earlier findings in literature (Zambrano et al, J Biol Chem. 1998 Aug 7;273(32):20128-33.; PMID: 9685356).

Positive Controls:

IMR32 neuronal cell lysate, HEK 293 cell lysate, or brain lysate.

Image: Western blot analysis using FE65 and HEK 293 cell lysates



FOR RESEARCH USE ONLY

NEUROMICS REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. v3-01/13

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail pschuster@neuromics1.com