

Patched tumor suppressor

Data Sheet

Catalog Number: RA19021

Product Type: Rabbit IgG

Immunogen Sequence: KKINYQVDSWE (a.a. 269-279) of the mouse PTCH protein.

Format: Affinity Purified liquid

Species Reactivity: Human, Mouse

Applications: Western blotting

Dilutions: WB 1:500-1:1000

Storage: Aliquot as undiluted antibody and immediately place at -20°C. Working aliquot may be stored at +4-8°C for up to 1 month. Liquid may have become trapped in top of vial during shipping. Centrifugation of vial is recommended before opening. Stable for at least 1 year at -20°C. Repeated freeze/thaw cycles compromise the integrity of the antiserum.