



---

**Catalog Number:** P16100

**Format:** liquid

**Description:** Metabotropic Glutamate Receptor 5  
(mGluR5) control peptide.

**Concentration:** 1 mg/ml

---

### Application Notes

Used for absorption control of antiserum binding. Antibody should be pre-incubated with peptide prior to incubation with tissue/sample. We suggest perincubation for 1 hour at room temperature. Dilution of 1:100 – 1:1000 is recommend to block antigen/antibody binding. However, optimal working concentration should be determined by the investigator using a range of peptide dilutions to eliminate or reduce immunoreactivity.

#### *FOR RESEARCH USE ONLY*

NEUROMICS ANTIBODIES REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.