



ORL-1 Peptide

Data Sheet

Catalog Number: P14133

Product Type: Peptide

Immunogen Sequence: CKTSETVPRPA

Molecular Weight: 1188 daltons

Description: Control Antigen for use in absorption control

Format: Liquid-100 ul.
Concentration of 1 mg/ml.

Storage: For maximum stability, the peptide at -20° C. Most peptides will be stable in solution for several days at 4° C. Avoid repeated freeze/thaws.

Application Notes

Used for absorption control of antiserum binding. Antibody should be pre-incubated with peptide prior to incubation with tissue/sample. We suggest pre-incubation for 1 hour at room temperature. A 10^{-6} or 10^{-5} molar (M) peptide concentration is suggested to block antigen/antibody binding. However, optimal working concentration should be determined by the investigator using a range of peptide dilutions to eliminate or reduce immunoreactivity.

FOR RESEARCH USE ONLY

NEUROMICS REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. v1-06/08

www.neuromics.com

Neuromics • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail pshuster@neuromics1.com