

NEUROMICS

GABA-B-R2



Datasheet

Catalog Number:	GT15009	Host:	Goat
Product Type:	Affinity Purified Goat IgG	Species Reactivity:	Rat
Immunogen Sequence:	Recombinant rat N-terminal extracellular domain of Gamma-Aminobutyric Acid Type B Receptor subunit 2	Format:	Lyophilized 0.2 µm filtered solution in phosphate-buffered saline (PBS) with 5% trehalose.
Applications:	Immunohistochemistry: 2-15 µg/mL Immunocytochemistry: 2-15 µg/mL Western Blot: 0.1 – 0.2 µg/mL ELISA: 0.5 - 1.0 µg/mL Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.		
Storage:	Antibodies are stable for six months from date of receipt when stored at -20° C to -70° C. The antibody can be stored at 2° - 8° C for 1 month without detectable loss of activity. Upon receipt, antibody should be aliquotted and stored frozen at -20° C to -70° C in a manual defrost freezer. <i>Avoid repeated freeze-thaw cycles</i>		

Application Notes

Reconstitution: Reconstitute with sterile PBS.

Immunohistochemistry: Antiserum was used on fixed tissue.

This antibody can be used with the appropriate secondary reagents to detect rat GABA-B-R2 on tissues dissected from animals that were fixed by vascular perfusion with 4% paraformaldehyde/PBS (pH 7.4) and followed by perfusion with a 10% sucrose solution in 0.1 M phosphate buffer (pH 7.2).

FOR RESEARCH USE ONLY

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

V. 11/03

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail pshuster@neuromics1.com