

Data Sheet

StrepMAB-Classic HRP

Cat. No.: 2-1509-001; 2-1509-002

Version: 13.0

Lot No.: 1509-

Revision Date: 25.11.2020

| | | |
|----------------------|--|---------|
| Description | Murine monoclonal antibody purified conjugated to horseradish peroxidase (HRP, 3-3.5 mols HRP/mol IgG). It is suitable for detection of Strep-tag®II and Twin-Strep-tag® proteins. | |
| Subclass | IgG1 | |
| Form | Suspension in PBS containing 5 mg/ml bovine serum albumin, approx. 0.05% ProClin®300 and stabilizing agents. | |
| Amount | 2-1509-001: 75 µl 2-1509-002: 20 µl | |
| Concentration | 1.5 mg/ml | |
| Application | ELISA/detection in coated microplates | 1:10000 |
| | Western blot, chromogenic reaction | 1:4000 |
| | Western blot, chemiluminescence reaction | 1:32000 |
| Stability | 24 months after shipping | |
| Storage | recommended: 2-8 °C | |
| Shipping | room temperature | |
| Hazards | Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided. | |



Please download all up-to-date manuals, protocols and other material from
<http://www.iba-lifesciences.com>.

For research use only

Trademark information

The owners of trademarks marked by “®” or “TM” are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.