



Catalog Number:	RA47000	Host:	Rabbit
Product Type:	Rabbit Monoclonal IgG	Species Reactivity:	Human
Immunogen Sequence:	Recombinant Human Clusterin/Apolipoprotein J	Format:	Liquid with PBS buffer, pH 7.4 with 0.02% Sodium Azide

Applications: Immunohistochemistry: 1:2-3

Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.

Storage: The product can be stored undiluted for several weeks at 4°C. Dilute only immediately before use. Aliquot and store at -20°C long term. Avoid freeze thaw-cycles.

Application Notes

Description/Data:

Clusterin/Apolipoprotein J Antibody detects for clusterin, a protein associated with the clearance of cellular debris and apoptosis which is encoded by the human gene, CLU. Used in a panel with other antibodies, detection of clusterin is useful for differentiating systemic anaplastic large cell lymphoma from classical Hodgkin's disease, or for the determination of follicular dendritic cell tumors.

FOR RESEARCH USE ONLY

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.-V2/08/2012

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail: pshuster@neuromics.com