




---

<b>Catalog Number:</b>	RA21026	<b>Host:</b>	Rabbit
<b>Product Type:</b>	Rabbit Polyclonal IgG	<b>Species Reactivity:</b>	Human, Mouse, and Rat
<b>Immunogen Sequence:</b>	Somatostatin conjugated to BSA.	<b>Format:</b>	Lyophilized from 10 mM PBS.
<b>Applications:</b>	Immunohistochemistry: 5-10 ug/ml Immunofluorescence: 5-15 ug/ml  Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.		
<b>Storage:</b>	The product can be stored as supplied for up to 12 months at 2°C-4°C. After reconstitution, aliquot and store at -20°C for higher stability or at 4°C with an appropriate antibacterial agent. Avoid freeze thaw-cycles.		

---

### Application Notes



#### Description/Data:

Somatostatin is secreted by neuroendocrine neurons in hypothalamus and interacts with somatostatin receptors inhibiting secretion of growth hormone.

*Images: Immunohistochemical staining of Somatostatin in rat brainstem medulla. 4% paraformaldehyde fixed rat brain 10 µm thick cryostat sections were incubated with rabbit polyclonal antibody to Somatostatin at 10 µg/mL overnight at 4°C followed by incubation with Donkey anti-rabbit Rhodamine Red (pseudo-coloured yellow) conjugated secondary antibodies at 1:200. Cell nuclei were counterstained with DAPI.*

**FOR RESEARCH USE ONLY**

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RSKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.-V2/08/2012

---

[www.neuromics.com](http://www.neuromics.com)

Neuromics Antibodies • 5325 West 74<sup>th</sup> Street, Suite 8 • Edina, MN 55439  
phone 866-350-1500 • fax 612-677-3976 • e-mail: [pshuster@neuromics.com](mailto:pshuster@neuromics.com)