



Catalog Number: PR27281

Product Type: Recombinant Protein

Source: E. Coli

Amino Acid Sequence: MGSSHHHHHH SGLVPRGSH MGSHMELLCH EVDPVRRRAVR DRNLLRDDR V LQNLITIEER YLPQCSYFKC VQKDIQPYMR RMVATWMLEV CEEQKCEEEV FPLAMNYLDR FLAGVPTPKS HLQLLGAVCM FLASKLKETS PLTAEKLCIY TDNSIKPQEL LEWELVVLGK LKWNLAAVTP HDFIEHILRK LPQQRKLSL IRKHAQTFIA LCATDFKFAM YPPSMIATGS VGAAICGLQQ DEEVSSLTCD ALTELLAKIT NTDVDCCLKAC QEQIEAVLLN SLQQYRQDQR DGSKSELD QASTPTDVRD IDL.

Description/Molecular Mass: CCND2 is a member of the highly conserved cyclin family, whose members are characterized by a dramatic periodicity in protein abundance through the cell cycle. Cyclins serve as regulators of CDK kinases. Various cyclins demonstrate distinct expression and degradation patterns which contribute to the temporal coordination of each mitotic event. CCND2 forms a complex with and acts as a regulatory subunit of CDK4 or CDK6, whose activity is vital for cell cycle G1/S transition. CCND2 interacts with and is involved in the phosphorylation of tumor suppressor protein Rb. High level expression of the CCND2 protein is observed in ovarian and testicular tumors. In addition, knockout studies of the homologous gene in mouse suggest the vital roles of the CCND2 gene in ovarian granulosa and germ cell proliferation.

CCND2 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 313 amino acids (1-289 a.a.) and having a molecular mass of 35.6kDa.

CCND2 is fused to a 24 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Purity: Greater than 85.0% as determined by:
(a) Analysis by SDS-PAGE.

Format: CCND2 protein solution (1mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 10% glycerol, 2mM DTT and 1M Urea.

Storage: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage, it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

FOR RESEARCH USE ONLY

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

03/08v1

www.neuromics.com

Neuromics • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail pshuster@neuromics.com