



Smooth Muscle Basal Medium

Catalog #: SMCM001

Size: 500ml

Storage: Store the medium at 2-8°C.

GENERAL INFORMATION

Smooth Muscle Basal Medium is a complete medium designed for the culture of smooth muscle cells. It was tested and optimized with smooth muscle cell growth and proliferation in vitro. It is formulated for use with 5% CO₂ and 95% air in a humidified incubator.

Smooth Muscle Basal Medium are subjected to comprehensive quality control tests using primary human cells. Approved in-house lots of Smooth Muscle Basal Medium are used as a reference. In addition, every lot of media is tested for the absence of microbial contaminants (fungi, bacteria, mycoplasma).

**Smooth Muscle Basal Media must be used for in vitro research purposes only. Our products are not authorized for human use, for in vitro diagnostic procedures, or for therapeutic procedures.*

Shipped in a Dry Ice Package.

FOR RESEARCH USE ONLY

NEUROMICS REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. v1-09809

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail pshuster@neuromics.com