

Data Sheet

Protein A Agarose

50% suspension

Cat. No.: 6-2010-001, 6-2010-002, 6-2010-005,
6-2010-025, 6-2010-100

Version: 3.1
Revision Date: 03.03.2020

Lot No.:

Description	Protein A Agarose resin for the purification of monoclonal- and polyclonal antibodies
Target substance	Polyclonal and monoclonal antibodies (IgG), bound via the Fc-region
Form	50 % suspension in buffer, pH 7.4: PBS pH 7.4, 1 mM EDTA, 0.01 % sodium azide.
Matrix	Rigid, highly cross-linked agarose
Agarose Average particle size ¹(Dv50)	45 µm
Ligand	Alkaline stable recombinant protein A expressed in <i>E. coli</i> using animal-free medium
Dynamic binding capacity ²(DBC10%)	≥70 mg hu IgG/mL resin at residence time 6 min ≥50 mg hu IgG/mL resin at residence time 2.4 min ≥25 mg hu IgG/mL resin at residence time 1 min
Recommended flow rate (20°C)	250 cm/h
Max flow rate	500 cm/h
Maximum back pressure	0.3 MPa, 3bar, 43 psi
Recommended pH working range	2 to 13 (normal working conditions)
Cleaning-In Place stability	No significant change in performance after 100 cycles with 0,5M NaOH (1cycle= 15 min contact time) use at room temperature
Chemical stability	Stable to commonly used aqueous buffers in Protein A chromatography
Storage	recommended: 2- 8°C
Shipping	room temperature
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.

¹Median particle size of the cumulative volume distribution.

²Dynamic binding capacity was determined with 1 mg/ml monoclonal hu IgG1 using a 1 ml column at 1, 2.4 and 6 minutes residence time.

For research use only

Trademark information

The owners of trademarks marked by "®" or "TM" are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.



Go digital and help the environment. Please download all up-to-date manuals, protocols and other material from <https://www.iba-lifesciences.com>.