

Data Sheet

pLSG-IBA45

Cat. No.: 5-4845-001

Version: 2.2

Lot No.: 4845-

Revision Date: 09.03.2020

Description	StarGate® Acceptor Vector designed for gene transfer into the polyhedrin gene locus of AcMNPV DNA by homologous recombination containing the following elements: <ul style="list-style-type: none"> • Polyhedrin promoter for high-level expression in insect cells. • Co-transfection with BacPAK6 linearized AcMNPV DNA (Clontech) or with circular <i>flashBAC</i> modified AcMNPV DNA (Oxford Expression Technologies) allows the generation of recombinant baculovirus at very high efficiency through reconstitution of an essential gene (ORF 1629) and elimination of wild-type virus to great extent. • Ampicillin resistance and ColE1 origin of replication (pUC) for propagation in <i>E. coli</i>. • The expressed recombinant protein will be localized in the cytoplasm.
Affinity tag	The recombinant protein will contain two affinity tags: <ol style="list-style-type: none"> 1. Strep-Tactin affinity tag (Strep-tag® II) for the purification of recombinant protein via Strep-Tactin resins. The Strep-tag® II is fused to the N-terminus of the recombinant protein. 2. 6xHistidine-tag for the purification of recombinant protein via Ni-NTA resins. The 6xHistidine-tag is fused to the C-terminus of the recombinant protein.
Resistance	Ampicillin
Form	5 µg, dissolved in 20 µl TE buffer, pH 8,0: 10 mM Tris-HCl, 1 mM EDTA
Concentration	250 ng/µl
Stability	12 months after shipping
Storage	recommended: 2-8 °C for frequent usage, -20 °C for long-term storage
Shipping	room temperature
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.

Note: The sequences have been compiled from information in the sequence database, published literature, and other sources, together with partial sequences obtained by IBA, however, the vectors have not been completely sequenced.



Go digital and help the environment. Please download all up-to-date manuals, protocols and other material from <http://www.iba-lifesciences.com>.

For research use only


Important licensing information

This product is based on StarGate, Strep-tag and 6xHistidine-tag technologies covered by intellectual property (IP) rights and on completion of the sale IBA grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at info@iba-lifesciences.com or at IBA GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.

Trademark information


The owners of trademarks marked by “®” or “TM” are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Expression cassette of pLSG-IBA45



LacP/Z cassette = contains LacZ alpha fragment under control of a separate promoter, which allows alpha complementation of *LacZ* mutations such as *LacZΔM15* as in *E. coli* DH5α or TOP10.

your protein = after StarGate cloning using *Esp3I* your gene of interest will be located here



Features	from bp	to bp	Sequencing primer
AcMNPVseq II	1	1395	LSG-Primer-for (Cat. No. 5-0000-161)
M13 ori	1447	1920	
Ampicillin resistance gene	2251	3111	5' - TAACCATCTCGCAAATAAATAAG -3'
ColEI ori	3259	3902	LSG-Primer-rev (Cat. No. 5-0000-162)
AcMNPVseq I	4211	5357	
Polyhedrin promoter	5286	5355	5' - CAACGCACAGAATCTAGCGC -3'
forward primer binding site	5286	5308	
<i>Strep-tag</i> II	5369	5410	
LacZ alpha fragment	5639	6040	
6xHistidine-tag	6104	6130	
reverse primer binding site	6196	6215	
total vector length		6215	