



IBA Lifesciences GmbH
Rudolf-Wissell-Str. 28
37079 Goettingen
Germany
Tel.: +49 (0) 551-5 06 72-0
E-mail: info@iba-lifesciences.com
www.iba-lifesciences.com

Data Sheet

MEXi-293E cells

Cat. No.: 2-6001-001

Version: 3.0
Revision Date: 02.06.2021

Description	MEXi-293E cells are human embryonic kidney (HEK) cells derived from the 293 cell line and used for propagation of plasmids and expression of recombinant proteins in mammalian cells. The cell line expresses the EBNA-1 protein for episomal replication of <i>oriP</i>-harboring plasmids. The cell line is adapted for growth in suspension and for use with the MEXi-Mammalian Expression System.
Cell type	Kidney (embryonic)
Growth properties	Suspension cells
Biosafety level	S1 (Germany), BSL-2 (USA) (Please check which level is applicable for your country).
Media adaptation	MEXi-CM and MEXi-TM medium
Form	cryopreserved
Amount	1.5 ml
Stability	6 months after shipping
Storage	store in cryo storage system at – 196 °C to -140 °C
Shipping	Dry ice
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.

For research use only

Trademark information

The owners of trademarks marked by “®” or “TM” are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale IBA Lifesciences grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at info@iba-lifesciences.com or at IBA Lifesciences GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.