

Data Sheet

MEXi-TM transfection medium

Cat. No.: 2-6011-010

Lot No.: 6011-

Version: 2.2

Revision Date: 27.02.2020

Description	MEXi-TM transfection medium is a medium especially developed for transient transfection and growth of MEXi-293E cells.
Medium Properties	Serum-free, chemical defined
Endotoxin (units)	<10 EU/ml
Media adaptation	is adapted to the MEXi-293E cell line
Form	liquid
Amount	1000 ml
Stability	6 months after shipping
Storage	recommended: 2- 8 °C, protect from light
Shipping	room temperature
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.



Go digital and help the environment. Please download all up-to-date manuals, protocols and other material from <http://www.iba-lifesciences.com>.

For research use only

Important licensing information

This product is based on StarGate technology covered by intellectual property (IP) rights and on completion of the sale IBA grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at info@iba-lifesciences.com or at IBA GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.

Trademark information

The owners of trademarks marked by "®" or "TM" are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.