

Data Sheet

10x Buffer BXT

Strep-Tactin®XT Elution Buffer with Biotin

Cat. No.: 2-1042-025

Lot No.: 1042-

Version: 2.1

Revision Date: 27.02.2020

Description	Buffer BXT is a Strep-Tactin®XT Elution Buffer with Biotin for competitive displacement of purified recombinant Strep-tag®II or Twin-Strep-tag® fusion proteins from immobilized Strep-Tactin®XT .
Components (10x concentrated solution)	Buffer BXT (Strep-Tactin®XT Elution Buffer with Biotin, 10x), pH 8.0 1 M Tris-Cl 1,5 M NaCl 10 mM EDTA 0,5 M Biotin
Stability	12 months after shipping
Storage	recommended: 2-8 °C
Shipping	room temperature
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. See provided Material Safety Data Sheet for details.
Working solution	Dilute 1 volume of concentrated solution (10x) with 9 volumes of demineralized or distilled water.



Go digital and help the environment. Please download all up-to-date manuals, protocols and other material from <http://www.iba-lifesciences.com>.

For research use only

Trademark information

The owners of trademarks marked by "®" or "TM" are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.